

PRODUCT SPECIFICATION
-----------------------

<b>Product name</b>	Male/Female/ Intermittent Injection Luer Locks
<b>Range</b>	Universal
<b>Classification</b>	Class IIa
<b>Marking</b>	Product name, Universal Logo, Size, Code Ref, Lot Number, Expiry date and CE Mark
<b>Material</b>	ABS plastic
<b>Packaging</b>	Pouch, dispenser, carton
<b>Sterilization</b>	EO Sterile

<b>Primary use</b>	Luer lock caps can be used for closing all standard male/female IV connections and syringes, but are not suitable for feeding syringes
--------------------	--

Product code	Colour	Quantity/inner box
UN877LF	Yellow	100
UN940	Red	100
UN940BL	Blue	100
UN940W	White	100
UN960	Yellow	100

<b>Quality certification</b>	ISO 13485
<b>Storage conditions</b>	Store in a cool dry place away from heat and direct sunlight
<b>Shelf life</b>	5 Years- if stored as recommended above
<b>Additional information</b>	Syringe caps provide a quick and convenient way to cap all IV syringes. Latex free

*Medical Device and Reporting System:  
According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

<b>REVISION:</b>	3.0	<b>DATE</b>	01/04/2013	<b>AUTHOR</b>	H Cocking
<b>NOTES:</b>	Updated to new format and periodic review				
<b>REVISION:</b>	4.0	<b>DATE</b>	5/04/2015	<b>AUTHOR</b>	F Li
<b>NOTES:</b>	Update the sterilization method				
<b>REVISION:</b>	5.0	<b>DATE</b>	02/06/2017	<b>AUTHOR</b>	R.Zhang
<b>NOTES:</b>	Periodic review, DSN format and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: <a href="http://www.shermond.com">www.shermond.com</a>	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--