

PRODUCT SPECIFICATION

<b>Product name</b>	Latex Examination Gloves
<b>Range</b>	Premier Protector
<b>Type</b>	Powder free, non-sterile
<b>Marking</b>	Product name, Premier Logo, Size, Code Ref., Lot Number, Expiry date, CE Mark, latex content and single use symbols
<b>Shape</b>	Standard cuff length and ambidextrous.
<b>Material</b>	Latex
<b>Colour</b>	Natural

Size	Length (mm) (according to EN455-2)	Palm width (mm) (according to EN455-2)	Product code
<b>XS</b>	240	76 ± 3	P7900
<b>S</b>	240	84 ± 3	P7901
<b>M</b>	240	94 ± 3	P7902
<b>L</b>	240	105 ± 3	P7903

<b>Primary use</b>	In the medical field to protect patient and user from cross-contamination when conducting medical examinations, diagnostic and therapeutic procedures and for handling contaminated medical material. Also, to reduce exposure between the wearer and external hazards such as chemicals and/or microorganisms.
--------------------	---

<b>Packaging</b>	Dispenser box (DB) XS; S; M; L.	100	items
	Outer Carton (OC)	10	DB

<b>Residual powder</b> (according to EN 455-3)	< 2 mg/glove
<b>Residual protein</b>	< 20ug/g
<b>Residual accelerator</b>	< 0.15%

<b>Thickness</b> (minimum at palm)	0.15 mm
---------------------------------------	---------

<b>Force at break</b> (according to EN 455-2)	≥ 6.0 N (during shelf life)
	≥ 6.0 N (after challenge testing)

<b>Quality certification</b>	ISO 13485
<b>Conformity to guidelines</b>	Class 1 - Medical Device Directive 93/42/EEC Cat. III - Personal Protective Equipment 89/686/EEC
<b>Conformity to standards</b>	EN 455 parts 1,2 & 3 BS EN 374-2:2003 micro-organism hazard BS EN374 -3:2003 chemical hazard Performance level 6 - Sodium Hydroxide 40% EN420:2003 +A1:2009
<b>Storage conditions</b>	Store in a cool dry place away from heat and direct sunlight
<b>Shelf life</b>	5 Years if stored as recommended above
<b>Additional information</b>	These single chlorinated gloves have an enhanced textured surface to provide excellent touch sensibility and grip, ideal when working with fluids, instruments or wet glassware.

*Medical Device and Reporting System:*

*According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

<b>REVISION:</b>	3.0	<b>DATE</b>	23/07/2015	<b>AUTHOR</b>	J Bertran
<b>NOTES:</b>	Periodic review, DSN format change and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: <a href="http://www.shermond.com">www.shermond.com</a>	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--