

DATA SHEET: DSN341

## PRODUCT SPECIFICATION

| Product name            | Latex Examination Gloves  |  |  |  |
|-------------------------|---|--|--|--|
| Range                   | Premier Prestige  |  |  |  |
| Туре                    | Powder free, sterile  |  |  |  |
| Marking                 | Product name, Premier Logo, Size, Code<br>Ref., Lot Number, Expiry date, CE Mark,<br>latex content, sterile and single use<br>symbols |  |  |  |
| Shape                   | Standard cuff length and ambidextrous.  |  |  |  |
| Colour                  | Natural   |  |  |  |
| Sterilization<br>Method | EO  |  |  |  |

| Size | Length<br>(mm)<br>(according to<br>EN455-2) | Palm width<br>(mm)<br>(according to<br>EN455-2) | Product code |
|------|---|---|--------------|
| S    | 240   | 84 ± 3  | P3301        |
| М    | 240   | 94 ± 3  | P3302        |
| L    | 240   | 105 ± 3   | P3303        |
|      |   |   |              |
|      |   |   |              |

In the medical field to protect patient and user from cross-contamination when conducting medical **Primary use** examinations, diagnostic and therapeutic procedures and for handling contaminated medical material

| Dackaging | Dispenser box (DB) XS; S; M; L. | 50 | pairs |
|-----------|---------------------------------|----|-------|
| Packaging | Outer Carton (OC)               | 4  | DB    |

This product complies with **Commission Directive** 2002/72/EC (contact with food)



| Residual powder         | < 2 mg/glove |  |
|-------------------------|--------------|--|
| (according to EN 455-3) | < 2 mg/giove |  |

| Thickness         | 0.15 mm |
|-------------------|---------|
| (minimum at palm) |         |

| Force at break (according to EN | ≥ 6.0 N (during shell life)       |  |  |
|---------------------------------|-----------------------------------|--|--|
| 455-2)                          | ≥ 6.0 N (after challenge testing) |  |  |

| Quality certification   | ISO 13485   |  |  |
|---|---|--|--|
| Conformity to guidelines  | Class 1s - Medical Device Directive 93/42/EEC                                     |  |  |
| Conformity to standards   | EN 455 parts 1,2 ,3 & 4   |  |  |
|   | BS EN374-2:2003 micro-organism hazard   |  |  |
|   | BS EN374-3:2003 chemical hazard Performance level 6 - Sodium Hydroxide 40%        |  |  |
| Storage conditions  | Store in a cool dry place away from heat and direct sunlight                      |  |  |
| Shelf life 5 Years if stored as recommended above   |   |  |  |
|   | Double chlorinated gloves with an enhanced textured surface to provide            |  |  |
| Additional information  | excellent touch sensibility and grip, ideal when working with fluids, instruments |  |  |
|   | or wet glassware.   |  |  |
| Medical Device and Reporting Sy   |   |  |  |
| According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices. |   |  |  |

| REVISON: | 3.0   | DATE | 23/07/2015 | AUTHOR | J Bertran |
|----------|---|------|------------|--------|-----------|
| NOTES:   | Periodic review, DSN format change and Shermond logo updated. |      |            |        |           |
| REVISON: | 4.0   | DATE | 09/11/2017 | AUTHOR | F Li      |
| NOTES:   | Updated sterilization method EO                               |      |            |        |           |

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: www.shermond.com Unit 6, Delta Park Industrial Estate Millmarsh Lane, London