

**DATA SHEET: DSN181** 

## PRODUCT SPECIFICATION

Product name	Softpore Surgical tape		
Range	Premier		
Marking	Each box is marked with Product name, Premier Logo, Size, Code Ref., Lot Number, Expiry date, CE Mark, latex free and single use symbols. Each roll has inner core marked with Premier Logo		
Colour	White		

Product code	Size		
Product code	Width (cm)	Length	
		(m)	
5080N	2.5	9.1 ± 5%	
5082N	5	9.1 ± 5%	
5084N	10	9.1 ± 5%	
5086N	15	9.1 ± 5%	

Primary use	For use in holding surgical dressings in place on patients. It may also be used in other medical areas as
	determined appropriate by trained medical staff.

Packaging				
Product code	Inner	Outer		
5080N	1	48		
5082N	1	24		
5084N	1	12		
5086N	1	8		

Material
Non woven polyester fibre + E>S Fibre: 45g/m <sup>2</sup> (±10%)
Glue coating: Synthetic rubber (St – Isoprene, Copolymer) Petroleum
Hydrocarbon Resin

Quality certification	ISO 13485		
Conformity to guidelines Class 1 – Medical Device Directive 93/42/EEC			
Conformity to standards	The material conforms to British Pharmacopeia specifications		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	3 Years if stored as recommended above		
Softpore tapes are self adhesive one-way stretch, made from 100% no polyester with an acrylic adhesive. Porous and breathable.  Ideal for securing swaps, dressings, cannulae and catheters. The product is latex free. The product is non sterile.			
Medical Device and Reporting Sy	stem: of MDD, incidents caused by this device must be reported immediately to our offices.		

REVISON:	5.0	DATE	05/02/2016	AUTHOR	J Bertran
NOTES:	Periodic review, DSN format and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond

Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ