

DATA SHEET: DSN454

PRODUCT SPECIFICATION

Product name	Gauze Swab					
Range	Premier					
Marking	Premier Logo, Size, Quantity, Code Ref, Lot No, Expiry Date, CE Marking, Latex free Symbol,					
	Non-Sterile Symbol					
Material	100% Cotton Gauze: Type 13 Light					
	Warp threads: 72/10cm +/- 2					
	Weft threads: 56/10cm +/- 2					
	14g/2 ² minimum mass					
Packaging	Paper Wrap Pack, Shipper Carton					
Sterilisation	Non-sterile					



Primary use	Intended to be used to cover or clean wound or abrasions and absorb their exudates	
	or to apply topical medications.	

Product	Colour	Dimension (cm)	Ply	Qty – Wrap Pack	Qty – Shipper Carton			
White								
1650	White	5 x 5	8	100x100	10,000			
1650A	White	5 x 5	8	120x100	12,000			
1665N	White	5 x 5	12	50x100	5,000			
1665A	White	5 x 5	12	80x100	8,000			
1655	White	7.5 x 7.5	8	100x100	10,000			
1655A	White	7.5 x 7.5	8	100x60	6,000			
1670N	White	7.5 x 7.5	12	50x100	5,000			
1660N	White	10 x 10	8	50x100	5,000			
1660A	White	10 x 10	8	100x60	6,000			
1675N	White	10 x 10	12	50x100	5,000			
1675A	White	10 x 10	12	100x40	4,000			
7011N	White	10 x 10	16	24x100	2,400			
7102	White	10 x 10	32	24x50	1,200			
Green								
7002A	Green	7.5 x 7.5	8	60x100	6,000			
7003N	Green	7.5 x 7.5	12	50x100	5,000			
7004N	Green	10 x 10	8	50x100	5,000			
7005A	Green	10 x 10	12	40x100	4,000			
7005N	Green	10 x 10	12	50x100	5,000			
7006A	Green	10 x 10	16	40x100	4,000			
Blue								
7008N	Blue	10x 10	8	50x100	5,000			
7009N	Blue	10 x 10	12	50x100	5,000			
7009A	Blue	10 x 10	12	40x100	4,000			
7010N	Blue	10 x 10	32	24x50	1,200			

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare T: +44 (0) 1530 278 111 W: www.365healthcare.com Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ



DATA SHEET: DSN454

Conformity to Regulations	European Medical Device Regulation 2017/745					
Device Classification	Class I					
Conformity to standards	100% Cotton Absorbent gauze in compliance with EN 14079 / BP 1988					
	Standards					
	ISO 13485:2016					
Storage conditions	Store in a cool dry place away from heat and direct sunlight					
Shelf life	5 Years- if stored as recommended above					
Additional information	Latex free					
Medical Device and Reporting System:						
Incidents caused by this device must be reported immediately to our offices.						

REVISON:	13.0	DATE	02/06/2020	AUTHOR	A Gidley
----------	------	------	------------	--------	----------

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.