

PRODUCT SPECIFICATION

Product name	Urine Drainage Bag
Range	Universal
Marking	Product name, Universal Logo, Code Ref, lot number, expiry date, size, CE mark, sterile (for UN40), single use and latex free symbols.
Material	Polyvinyl chloride (PVC)
Volume	2 litres

Product code	Product type	Quantity per outer
UN40	With taps, sterile	50
UN2	Non sterile	2 x 250

Primary use	To collect urine. The bags are suitable for use with a urine drainage bag stand or hanger
--------------------	---

Product dimensions		
	UN40	UN2
Length of the PVC tube	950 ±10mm	950 ±10mm
Distance between outer holes	155 ±1mm	160 ±1mm
Distance between inner holes	75 ±1mm	75 ±1mm
Intern diameter connector	5.85 ±0.15mm	5.7 ±0.5mm
Intern diameter of Inlet tube	6.0 ± 0.15mm	5.0 ±0.5mm

Quality certification	ISO 13485
Conformity to guidelines	UN40 - Class Is – Medical Device Directive 93/42/EEC UN2 - Class I – Medical Device Directive 93/42/EEC
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	5 Years if stored as recommended above
Additional information	Closed system non drainable urine bags with universal connector. The product is leak, impact and pressure tested according to the manufacturer standards. UN2 - non-return valve and tear-open-to-empty. UN40 - sample port and lever tap. The products are latex free and single use.
<i>Medical Device and Reporting System: According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.</i>	

REVISION:	3.0	DATE	25/02/2016	AUTHOR	J Bertran
NOTES:	Update product specs, DSN format and Shermond logo.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: www.shermond.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--